



### Product Highlights

- Deluxe Procedure Kit
- Available in Both Pull and Push Techniques
- Universal-Sized Snare
- Multiple Sizes and Accessories
- Indicated for Pediatric and Adult Use

### PRODUCT FEATURES

#### High-Grade Silicone Dome-Bolstered PEG Tube

- Designed for performance and consistent flow
- May facilitate traction removal

#### Universal 1.9mm O.D. x 240cm Retrieval Snare

- Traverses scopes with a 2.0mm or larger working channel
- Intended to eliminate the cost of a specialty snare

#### Double-Tier Tray Packaging

- Provides kit organization and may promote procedural efficiency

#### External Bolsters

- Two external bolsters provided to accommodate physician preference and patient needs: Round Bolster and Roller Bar Bolster

#### Multiple Sizes and Accessories

- Multiple French sizes for broad patient coverage
- Separate medication and feeding ports
- C-Clamp designed to help reduce risk of tube reflux and leaking

### KIT CONTENTS

- |                                      |  |
|--------------------------------------|--|
| (1) Fenestrated Drape                | (1) Scissors                             |
| (1) No. 11 Blade Surgical Scalpel    | (1) Curved Hemostat                      |
| (1) Retrieval Snare                  | (2) External Bolsters                    |
| (1) Insertion Wire                   | (1) C-Clamp                              |
| (1) Guidewire                        | (1) Y-Port                               |
| (1) High-Grade Silicone PEG Tube     | (4) 2" x 2" (5cm x 5cm) Drainage Sponges |
| (4) 4" x 4" (10cm x 10cm) Gauze Pads | (1) Angiocath™ and/or Seldinger Needle   |



### Ordering Information

#### STANDARD PEG KITS

Order Number	Description	Unit
M00568381	20Fr (6.67mm) Standard PEG Kit, Pull	Box 2
M00568401	20Fr (6.67mm) Standard PEG Kit, Push	Box 2
M00568391	24Fr (8mm) Standard PEG Kit, Pull	Box 2
M00568411	24Fr (8mm) Standard PEG Kit, Push	Box 2

#### PEG ACCESSORIES

Order Number	Description	Unit
M00580641	20Fr (6.67mm) Y-Ports	Box 10
M00580651	24Fr (8mm) Y-Ports	Box 10
M00580731	20Fr (6.67mm) Round Bolsters	Box 10
M00580781	20Fr (6.67mm) and 24Fr (8mm) C-Clamps	Box 10

**European Distribution Centre – The Netherlands**  
T: +31 45 54 67 700 F: +31 45 54 67 800

**European Headquarters – Paris**  
T: +33 1 57 66 80 00 F: +33 1 57 66 84 99

**Argentina (Freephone)**  
T: 0800 555 2678 F: +5411 4896 8517

**Australia/New Zealand**  
T: +61 2 8063 8100 F: +61 2 9330 1404

**Austria**  
T: +43 1 60 810 F: +43 1 60 810 60

**Belgium (Freephone)**  
T: 0800 94 494 F: 0800 93 343

**Brazil**  
T: +55 11 5502 8500 F: +55 11 5103 2212

**Canada**  
T: +1 888 359 9691 F: +1 888 575 7396

**Chile**  
T: +562 445 4904 F: +562 445 4915

**China – Beijing**  
T: +86 20 8767 9791 F: +86 20 8767 9791

**China – Guangzhou**  
T: +86 20 8767 9791 F: +86 20 8767 9789

**China – Shanghai**  
T: +86 21 6141 5959 F: +86 21 6141 5900

**Colombia**  
T: +57 1 629 5045 F: +57 1 612 4761

**Czech Republic**  
T: +420 296 331 901 F: +420 296 331 935

**Denmark (Freephone)**  
T: 80 30 80 02 F: 80 30 80 05

**Eire (Freephone)**  
T: 1800 882 969 F: 1800 882 968

**Finland**  
T: +358 20 762 8882 F: +358 20 762 8883

**France**  
T: +33 1 39 30 49 00 F: +33 1 39 30 49 01

**Germany (Freephone)**  
T: 08000 723300 F: 08000 723319

**Greece**  
T: +30 210 9542 300 F: +30 210 9542 310

**Hong Kong**  
T: +852 2960 7100 F: +852 2563 5276

**Hungary**  
T: +36 1 456 30 40 F: +36 1 456 30 41

**India – Bangalore**  
T: +91 80 2212 4928/9 F: +91 80 2207 5153

**India – Chennai**  
T: +91 44 2220 1879 F: +91 44 2220 0716

**India – Delhi**  
T: +91 11 4243 2222 F: +91 11 2610 0808

**India – Mumbai**  
T: +91 22 4030 9165 F: +91 22 4040 9199

**Italy**  
T: +39 010 60 60 1 F: +39 010 60 60 200

**Korea**  
T: +82 2 3476 2121 F: +82 2 3476 1776

**Mexico**  
T: +52 55 5687 63 90 F: +52 55 5687 62 28

**Middle East/Gulf/North Africa**  
T: +961 1 805 410 F: +961 1 805 445

**The Netherlands**  
T: +31 30 602 55 44 F: +31 30 602 55 05

**Norway (Freephone)**  
T: 800 104 04 F: 800 101 90

**Poland**  
T: +48 22 435 14 14 F: +48 22 435 14 10

**Portugal**  
T: +351 21 381 25 40 F: +351 21 381 25 58

**South Africa**  
T: +27 11 840 8600 F: +27 11 463 6077

**South East Asia – Malaysia**  
T: +60 3 2283 3813 F: +60 3 2284 3813

**South East Asia – Philippines**  
T: +63 2 687 6994 F: +63 2 687 3047

**South East Asia – Singapore**  
T: +65 6418 8888 F: +65 6418 8899

**South East Asia – Thailand**  
T: +66 2 6380 100 F: +66 2 6380 400

**South East Europe**  
T: +30 210 9542 300 F: +30 210 9542 310

**Spain – Barcelona**  
T: +34 93 444 72 00 F: +34 93 405 90 45

**Spain – Madrid**  
T: +34 90 11 12 15 F: +34 91 319 50 03

**Sweden**  
T: +46 42 25 69 00 F: +46 42 25 69 69

**Switzerland**  
T: +41 32 626 57 00 F: +41 32 626 57 01

**Taiwan**  
T: +886 2 2747 7278 F: +886 2 2747 7270

**Turkey – Istanbul**  
T: +90 216 464 36 66 F: +90 216 464 36 67

**UK (Freephone)**  
T: 0844 800 4512 F: 0844 800 4513

**Uruguay**  
T: +59 82 900 6212 F: +59 82 900 6212

**Venezuela**  
T: +58 212 959 6275 F: +58 212 959 5328

# Boston Scientific

Delivering what's next.™

www.bostonscientific.com  
www.bostonscientific-international.com